

SHARPER IMAGE®

SHIATSU KNEADING AND ROLLING FOOT MASSAGER

Item No. 205365

User Guide



Thank you for purchasing the Sharper Image Shiatsu Kneading and Rolling Foot Massager. This device uses modern technology to replicate ancient Chinese shiatsu massage techniques to stimulate the vital points on your feet.

FEATURES

- 4 Preset massage modes
- 3 Auto target modes for toe / arch / sole
- Manual mode with variable speed and direction
- Optional, removable foot comfort pads for a less intense massage
- 15-minute automatic shutoff

BASIC OPERATION

- Plug the power cord directly into a standard AC outlet.
- Press the ON/OFF button to turn the unit on.
- Set the unit to your choice of massage programs.

USING THE CONTROL PAD

Auto Mode: Press the “Auto” button to cycle through four preset massage programs. Each program contains a specific massage sequence designed to stimulate the different pressure points on your feet, calves, heels and legs. [You will hear a beep each time a new massage program is selected.]

P1: Clockwise medium speed 3 min. / clockwise high speed 4 min. / counterclockwise medium speed 4 min. / counterclockwise high speed 4 min.

P2: Clockwise medium speed 2 min. / clockwise low speed 3 min. / counterclockwise medium speed 3 min. / counterclockwise low speed 2 min. / clockwise low speed 3 min. / clockwise high speed 2 min.

P3: Clockwise medium speed 3 min. / counterclockwise high speed 3 min. / clockwise low speed 3 min. / counterclockwise medium speed 3 min. / clockwise high speed 3 min.


P4: Clockwise high speed 4 min. / clockwise low speed 3 min. / counterclockwise medium speed 2 min. / counterclockwise high speed 3 min. / counterclockwise low speed 3 min.

Custom Mode: Press the “Custom” button to automatically target different areas of the foot [Toe, Arch and Sole].

Manual Mode: Press the “Manual” button to begin a personalized massage experience. Simply use the Speed and Direction arrows to find the perfect massage for your needs.

USING THE INCLUDED FOOT COMFORT PADS (OPTIONAL)

If you find the massage nodes are too intense for your feet, you may use the included Foot Comfort Pads. Simply place the pads in the massage ports [where the nodes meet your feet] for an additional layer of padding and comfort.



AUTOMATIC SHUTOFF

Please note: there is an automatic shutoff timer that turns off all massage programs after 15 minutes.

SAFETY WARNINGS

WARNING: To reduce the risk of burns, fire, electric shock or injury to people:

- Only use this appliance for its intended household use as described in this manual.
- Before using this device, inspect the unit and power cord to determine if there is any damage. If there is any damage, do not use this device.
- Do not use this device in or near water. Do not immerse in water.
- Not for use by children. Keep this device out of reach of children.

WARRANTY / CUSTOMER SERVICE

Sharper Image branded items purchased from SharperImage.com include a 1-year limited replacement warranty. If you have any questions not covered in this guide, please call our Customer Service department at 1 (877) 210-3449. Customer Service agents are available Monday through Friday, 9:00 a.m. to 6:00 p.m. ET.

SHARPER IMAGE®

The SHARPER IMAGE® name and logo are registered trademarks.

©Sharper Image. All Rights Reserved.